

# The 27-gene IO score is associated with molecular features and response to immune checkpoint inhibitors (ICIs) in patients with gastric cancer



Matthew G. Varga<sup>1\*</sup>, Catherine Cronister<sup>1</sup>, Tyler Nielsen<sup>1</sup>, Kim McGregor<sup>1</sup>, Doug Ross<sup>1</sup>, David R Hout<sup>1</sup>, Robert Seitz<sup>1</sup>, David Gandara<sup>2</sup> and Brock L. Schweitzer<sup>1</sup>

<sup>1</sup> Oncocyte Corporation, Nashville, TN USA <sup>2</sup> UC Davis Comprehensive Cancer Center, Sacramento, CA, USA. \* Correspondence: mvarga@oncocyte.com

## Introduction

- Gastric cancer (GC) is the 3rd leading cause of cancer-related death worldwide.
- A variety of biomarkers have been used to identify patients most likely to benefit from ICI therapies such as high PD-L1 expression, MSI-high, Epstein Barr virus (EBV) positive, or TMB high.
- ICIs + chemo have recently been approved as a first line treatment in advanced GC for patients with CPS > 5%.
- Among GCs, there is extensive molecular heterogeneity, tumor immune microenvironments (TIME), and patient factors that may be associated with outcomes
- Despite these potential biomarkers, most patients with advanced GC do not respond to ICI treatment.

## The 27-Gene IO Assay

- The CLIA-validated 27-gene immuno-oncology (IO) assay was developed to assess both inflammatory effector cell and surrounding cancer associated fibroblasts phenotypes in the TIME (Figure 1).




Figure 1: The IO score is a measure of both the quiescent and active states of the TIME




Table 1: Cohort Characteristics of TCGA and ACRG

Feature	TCGA (STAD)			ACRG		
	N	IO+ (n,%)	IO- (n,%)	N	IO+ n (%)	IO- n (%)
All Patients	135	59 (44%)	76 (56%)	294	125 (43%)	169 (57%)
Age			0.49			0.56
<64 years	55	22 (40%)	33 (60%)	151	67 (44%)	84 (56%)
≥ 64 years	80	37 (46%)	43 (54%)	143	58 (41%)	85 (59%)
Sex			0.097			0.90
Female	47	25 (53%)	22 (47%)	99	43 (43%)	56 (57%)
Male	88	34 (39%)	54 (61%)	195	82 (42%)	113 (58%)
Stage			0.40			0.05
I	18	10 (56%)	8 (44%)	30	19 (63%)	11 (37%)
II	59	25 (42%)	34 (58%)	92	40 (43%)	52 (57%)
III	55	24 (44%)	31 (56%)	94	38 (40%)	56 (60%)
IV	3	0 (0%)	3 (100%)	76	28 (37%)	48 (63%)
Lauren Type			0.26			1
Intestinal	97	45 (46%)	52 (54%)	141	57 (40%)	84 (60%)
Diffuse	23	7 (30%)	16 (70%)	133	58 (44%)	75 (56%)
Mixed	11	3 (27%)	8 (73%)	17	7 (41%)	10 (59%)
Tumor Location			0.89			0.24
Antrum	52	21 (40%)	31 (60%)	152	69 (45%)	83 (55%)
Body	39	19 (49%)	20 (51%)	104	37 (36%)	67 (64%)
Cardia	25	11 (44%)	14 (56%)	32	17 (53%)	15 (47%)
Fundus	17	7 (41%)	10 (59%)	6	2 (33%)	4 (66%)

Table 2: Molecular Features of TCGA and ACRG Cohorts

Feature	TCGA (STAD)			ACRG		
	N	IO+ n (%)	IO- n (%)	N	IO+ n (%)	IO- n (%)
EBV			p<0.001			p<0.001
EBV+	14	14 (100%)	0 (0%)	18	18 (100%)	0 (0%)
	121	45 (37%)	76 (63%)	251	94 (37%)	157 (63%)
MSI			p=0.02			p=0.016
MSI High	82	29 (35%)	53 (64%)	111	58 (53%)	53 (48%)
MSS	53	30 (57%)	23 (43%)	183	69 (38%)	114 (62%)
TMB			p<0.001			Data Not Available in ACRG
High (>10mut/MB)	36	26 (72%)	10 (28%)			
Low	99	33 (33%)	66 (67%)			
PDL1 mRNA Expression			p<0.001			p<0.001
High (>75th quartile)	34	28 (82%)	6 (18%)	74	61 (82%)	13 (18%)
Low	101	31 (31%)	70 (69%)	74	3 (4%)	71 (96%)

Table 3: Molecular features associated with IO score in a clinical cohort

Feature	N	IO+ n (%)	IO- n (%)
All Patients	59	25 (42%)	34 (58%)
EBV			p=0.15
Positive	5	4 (80%)	1 (20%)
Negative	40	15 (37%)	2 (63%)
Missing	14	6 (43%)	8 (57%)
MSI			p=0.79
MSS	39	16 (41%)	23 (59%)
MSI-High	20	9 (45%)	11 (55%)
PDL-1 (CPS%)			p=0.02
Low (<1%)	20	3 (15%)	17 (85%)
High	32	16 (50%)	16 (50%)
Missing	7	6 (85%)	1 (15%)
TCGA Classifier			p=0.04
EBV	4	4 (100%)	0 (0%)
GS	20	8 (40%)	12 (60%)
CIN	16	4 (25%)	12 (75%)
MSI-H	5	3 (60%)	2 (40%)
Missing	14	6 (42%)	8 (57%)
Best Response			p=0.04
CR	4	3 (75%)	1 (25%)
PR	15	10 (67%)	5 (33%)
SD	21	7 (33%)	14 (67%)
PD	19	5 (26%)	14 (74%)

Table 4: Response status by molecular feature

Feature	Responder CR/PR	Non-Responder PD/SD
All Patients	19 (32%)	40 (68%)
MSI type		
MSI	21 (60%)	8 (40%)
MSS	7 (18%)	32 (82%)
IO Score		
Positive	13 (52%)	12 (48%)
Negative	6 (18%)	28 (82%)
EBV		
Positive	5 (100%)	0 (0%)
Negative	7 (18%)	33 (83%)
Missing	7 (50%)	7 (50%)
PD-L1 1%		
≥ 1%	16 (50%)	16 (50%)
< 1%	3 (11%)	24 (89%)
TCGA Classifier		
EBV	4 (100%)	0 (0%)
GS	2 (10%)	18 (90%)
CIN	1 (6%)	15 (94%)
MSI-H	5 (100%)	0 (0%)
Missing	7 (50%)	7 (50%)




Figure 2: A) IO continuous score plotted by TCGA classifier status B) IO continuous score by objective response. Dotted line at 0.09, the binary threshold for IO positivity.

## Aim

To demonstrate that the 27 gene IO score, a TIME classifier, is associated with the existing molecular markers of gastric cancer and with objective response to ICI therapy in a clinical cohort.

## Methods

- RNA-seq expression data was obtained from 3 independent cohorts including TCGA (STAD), ACRG (GSE84437, GSE84426), and a clinical cohort with ICI response data (PRJEB25780, PRJEB40416)
- The 27 gene IO algorithm was applied to all available patient data to derive IO scores
- Fisher's exact test was used to examine the associations between IO score and clinical features and molecular subtypes.

## Results

- In TCGA, IO score was associated with the molecular features EBV, MSI, TMB, and PD-L1
- In ACRG, the IO score was significantly associated with EBV, MSI, and PD-L1
- IO positive found 36 of 37 EBV+ patients across all cohorts
- In the clinical cohort of 59 patients, IO score was associated with ICI response (p<0.05)
- Association between IO score and ICI response OR= 5.3 (95% CI: 1.3 to 23.92, p=0.01).

## Conclusions

Further studies are warranted to demonstrate that the 27-gene IO score may be a more comprehensive biomarker for assessing the TIME and provide complementary data to tumor-specific biomarkers for gastric cancers.